# LevoCit

#50300151

30.09.2014

Instructions

Anleitung

Mode d'emploi

# English

## Description

LevoCit is an off-white, two component acrylic mounting system with filler, consisting of one powder and one liquid. LevoCit is used when good edge retention and good workability properties are required.

LevoCit cures at room temperature.

## Material Safety Data

### *LevoCit Powder*

Please read the following carefully:

Contains: Dibenzoyl peroxide.  
May cause an allergic skin reaction.

### *LevoCit Liquid*

DANGER

Contains: Methyl methacrylate

H225: Highly flammable liquid and vapour.

H315: Causes skin irritation.

H317: May cause an allergic skin reaction.

H335: May cause respiratory irritation.

P210: Keep away from heat/sparks/open flames/hot surfaces. – No smoking.

P261: Avoid breathing dust/fume/gas/mist/vapours/spray.

P280h: Wear protective clothing and gloves.

P302+352: IF ON SKIN: Wash with soap and water.

P333+313: If skin irritation or rash occurs: Get medical advice/attention.

## Mounting of Specimens

### *Measuring of Components*

The ratio for LevoCit Powder:LevoCit Liquid is **20:10** by weight and **20:10** by volume. Weighing is recommended to obtain the best mounts and for reproducibility.

#### ***Important!***

Do NOT pour excess liquid back into the bottle as this will significantly reduce the shelf-life.

### *Mixing of Components*

- Measure out the amount of LevoCit Liquid into a mixing cup (max. 30 g / 30 ml).
- Measure out the amount of LevoCit Powder into a mixing cup (max. 60 g / 60 ml).
- Pour the powder into the liquid.
- Stir thoroughly for **45 s.**, until the mixture is homogeneous and without lumps.  
Mix carefully to avoid trapping air in the mixture.

### *Selecting the Mounting Cup*

All Struers mounting cups can be used.  
If possible, choose a mounting cup where the distance between the specimen and cup wall is approx. 5 mm / 0.2".

### *Pre-treatment*

To promote the best possible adhesion of LevoCit to the specimen ensure that the specimen is thoroughly cleaned and dry prior to mounting. Each surface must be free from grease, moisture and other contaminants.

### *Mounting*

- Place the cleaned specimen in the centre of the mounting cup.
- Mix the required amount of LevoCit.
- Pour the mixture carefully over the specimen; run the mixture slowly along the stirrer or the mixing cup handle.
- Ensure that no air remains trapped within cracks or cavities in the specimen by moving the specimen or by tapping the mounting cup on the table.

### **Curing**

LevoCit is a room temperature curing system. For a 30 mm/ 1¼" dia. mount, the curing time is approx. 20 minutes. We recommend keeping the mount in the mounting cup until the temperature is below 40°C/104°F.

### *After Curing*

After removal from the mounting cup, the mount is immediately ready for preparation.

# Technical Data

## Mixing ratio

	Components	Ratio by weight	Ratio by volume
LevoCit	LevoCit Powder	20	20
	LevoCit Liquid	10	10

**NB** For the best mount use weighing method.

## Density

LevoCit Powder: 0.95 g/ml

LevoCit Liquid: 0.97 g/ml

## Soluble in

LevoCit Liquid is soluble in Acetone

## Filler

The filling material in LevoCit is calcium carbonate.

## Curing

Can be prepared (ground and polished) after:	20 min.
Peak temperature *:	75°C / 167 °F
Pot life after mixing:	3½ min.

\*30 mm dia. mounting cup without specimen at 21°C / 70°F.

## Properties after Curing

Hardness: 84 Shore D

Density: 1.4 g/ml

Chemical resistance: LevoCit is resistant to the most common etching materials but can be coloured.

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**Struers A/S**